RF4

FOCUS DESIGN AND CONCEPT


RF-4KC1

Industrial camera manual

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Directions for use Catalogue

I . Outline structure drawing

II . The equipment parameters

The parameter types	The parameter value
Sensor Target surface size	1/2.8"
The sensor resolution	4k/3840*2160@30fps
Output interface	HDMI 4k@30fps
Mode of operation	button
Working voltage	DC 12V
Storage	Usb disk, capacity<=64GB
The lens type	Supports C/CS interfaces
Photograph	8 million, JPG
Video	3840*2160 avi
Cross curve	8 * 8 groups
File Manager	With
Parameter type	specific parameter

Tip:

Hold down the '1080 key' for 10 seconds to 1080P resolution

Hold down the '4K key' for 10 seconds to switch to 4K resolution

III. Packing list

1. One host
2. One DC 12V/ 1A power adapter

Note: if the above packing parts are missing, please contact the dealer or our company to ensure your rights and interests.

IV. Install and use

1. The microscope is connected to the host computer through C/CS interface, and tightened to prevent falling;
2. HDMI connection to monitor. Display should support 3840*2160/1920*1080 resolution;
3. Use USB interface to connect computer for UVC camera;
4. Connect the power adapter to the host through the DC port.
5. The TF card can be used with related functions: taking photos, video recording, playback, etc. If necessary, please insert it;
6. Switch on the power adapter, and the host indicator shows green for 10 seconds.

FAQ:

1. Problem: The power indicator does not change after power on.
Solution: Ensure that the power adapter is powered on properly and the DC connector is not oxidized.
2. Problem: The display does not respond after powering on.
Solution: Check the following: 1) The HDMI interface is properly connected to the display interface. 2) Display power, switch, input interface selection is normal; 3) Display supports 3840*2160/1920*1080P resolution;

V. Operation instructions

5.1 Key Operation Instructions

- 1) Button 1: When there is a storage device, short press to capture, long press to record and stop recording;
- 2) Press 2, the device will switch in the following states: menu, cross line, playback (when the storage device is plugged in);
- 3) When the device is in different states, the functions of buttons 3, 4, 5 and 6 are as follows:
 - 1))Menu status: 3 and 4 are used to adjust the highlighted button up and down, 5 and 6 are used to perform key functions or adjust function parameters;
 - 2))Crosshead status: 3 and 4 are used to adjust the highlighted crosshead. 5 and 6 are used to adjust the position of the crosshead.
 - 3))Playback state: 3,4 are used to flip up and down in single screen or multi-screen; 5 are used to switch between single screen and multi-screen; 6 are used to display a single screen and select a picture or video in multi-screen state;

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The home page	
Image parameters	Adjust the image effect;
Function setting	Adjust language, cross line and other functional parameters;
Device management	Set usb device and resolution.;
Pop up the disk	Eject storage media such as usb flash drive;
Exit	Exit menu;

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Image parameters	
Acutance	Adjust the sharpness of the image;
luminance	Adjust the brightness of the image;
Contrast	Adjust the contrast of the image;
Tonal	Adjust the warm and cold colors of the image;
saturation	Adjust image color saturation;
Manual exposure	Switch the device exposure to manual exposure or automatic exposure;
Auto white balance	Switch the white balance of the device to manual or automatic white balance;
One key white balance	In manual white balance mode, the lens locks the white balance parameters with reference to the white sample;
Red	In manual white balance mode, red gain value;
Green	In manual white balance mode, green gain value;
Blue	In manual white balance mode, blue gain value;
Return	Return to main menu;

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Function Setting	
Language: Chinese	Set the device language to Chinese
Language: English	Set the device language to English
rotating	Rotate the image 180 degrees
Current Line group	Sets the current line group
Number of vertical lines	Sets the number of vertical lines in the line set currently in use
Number of horizontal lines	Sets the number of horizontal lines in the line group currently in use
Line width	Set the width of the crosshead
Button left and right to adjust the color of the line	Set the device in crosshair mode, The left and right buttons (5, 6) function to adjust the color
Adjust line position by pressing left and right buttons	Set the device in crosshair mode, Left and right (5, 6) buttons function to adjust line position
Return	Returns to the main menu

5.5 Sub-menu (Device Management)


设备管理	
3840X2160@30FPS	[NO]
2560X1440@30FPS	[NO]
1920X1080@60FPS	[YES]
1280X1024@60FPS	[NO]
1280X720@60FPS	[NO]
1024X768@60FPS	[NO]
恢复出厂设置
升级
返回
Version:15:23:56, Dec 4 2021	
菜单模式	

Equipment management	
Resolution	3840x2160; 2560X1440; 1920x1080 ; 1280x1024; 1280x720; 1024x768;
Restore factory Settings	The system restarts after factory Settings are restored;
Upgrade	Upgrade can be performed based on the upgrade package provided by the vendor
Return	Return to the main menu
Current Version information	

VI. Precautions for Customer use

1. Insert a usb flash drive before taking photos or videos.
2. Disconnect the power supply when the device is not in use to prolong the service life of the device.
3. Do not insert the power adapter whose voltage exceeds the device limit.
4. Avoid violent collision.

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